

# FastClamp™

## Endoscopic Clamping System

FastClamp™ is a high specification, cost-effective endoscopic clamping solution. It ensures quick, stable, unobtrusive instrument and scope fixation.

The device is comprised of an instrument clamping arm and table clamp.


### FastClamp™ clamping arm:

- + Has multiple degrees of freedom yet can be securely locked to allow instruments to be positioned accurately and quickly
- + Securely grips either 5mm or 10mm instruments by way of the 'soft jaws' which also minimise the risk of instrument damage
- + Uniquely designed to allow one or two arms to be used with each table clamp
- + Two arm lengths available: standard and bariatric, ideally complementing PretzelFlex™

For further detail about FastClamp™ please contact:

[sales@surginno.co.uk](mailto:sales@surginno.co.uk)

+44(0)113 230 7597



### FastClamp™ table clamp:

- + Easily positioned with a wide range of options of adjustment for secure yet rapid 'over the drape' clamping to any standard table rail
- + The clamp handle's toggle action ensures a rigid connection to the operating table
- + Whole system can be released in seconds



# FastClamp™

## Endoscopic Clamping System

### FastClamp™ clamping arm options

Product Name	Product Code
FastClamp™ System	A-1000
Extended FastClamp™ System	A-1001

### FastClamp™ table clamp options

Product Name	Product Code
Table Clamp	TC-1000-01
Instrument Clamp	IC-1000-01
Extended Instrument Clamp	IC-1001-01

### Bar options

Product Name	Product Code
300mm Bar	BC-1030-01
500mm Bar	BC-1050-01
700mm Bar	BC-1070-01
Bariatric Bar	BC-1090-01



creative solutions for minimally invasive surgery

[www.surginno.com](http://www.surginno.com)

Leaflet Part Number: 103-9001 Issue1  
Change Request Number: CR3571