


# YelloPort Elite®

Resposable® Laparoscopic Port Access System


RESPOSABLE®  
PRODUCT

## Universal Seal

- Accommodates 5-12mm instruments
- Resposable® - **The cost effective port access system**
- Lightweight and low profile design



## Single-use Universal Seal





## Product Codes

Product Name	Product Code	Diameter
Universal Seal	EA512US	5-12mm

	Product Name	Product Code	Working Length
10mm	Cannula	EC1010520	105mm
	Pencil Point Trocar	ET1010503	105mm
	Blunt Trocar	ET1010504	105mm
	Cannula	EC1015020	150mm
	Pencil Point Trocar	ET1015003	150mm
	Blunt Trocar	ET1015004	150mm

	Product Name	Product Code	Working Length
12mm	Cannula	EC1210520	105mm
	Pencil Point Trocar	ET1210503	105mm
	Blunt Trocar	ET1210504	105mm
	Cannula	EC1215020	150mm
	Pencil Point Trocar	ET1215003	150mm
	Blunt Trocar	ET1215004	150mm

For further details about YelloPort Elite® please contact:

[sales@surginno.co.uk](mailto:sales@surginno.co.uk)

+44(0)113 230 7597

Brochure Part Number: 479-9001 Issue: 2  
Change Request: CR4551